

Balanced Flow Multi-Station Manifolds

THE ONLY multi-station delivery system that distributes equal rates of fresh anesthetic gas to each port independent of what accessory is used.



Features

- Designed for multiple mouse, rat and other small rodent procedures.
- Can accommodate up to six anesthesia stations/accessories.
- Equal fresh gas flow to each port independent of accessory used.
- Includes 15mm common outlets for easy connection to anesthesia accessories.
- Mounts onto any Scivena Scientific Anesthesia Machine and most other manufacturer's anesthesia machines with available mounting brackets.



Available Configurations:

RES535- 2 Port
RES536- 4 Port
RES537- 6 Port