


Multi-Station Rodent Procedure Board

SCIVENA
SCIENTIFIC
ANESTHESIA INNOVATIONS

RES3150

The only multiple rodent procedure system that assures equal fresh gas distribution to each subject.



****Shown with M3000 Anesthesia Machine, optional vaporizer, and Balanced Flow Multi-Station Manifold**

Features:

- Six numbered subject positions- plugs included for un-used positions
- 22mm Evacuation post for passive waste gas management
- Adjustable tie-downs (extra material included)
- Includes six industry standard RES500 coaxial Rat and Mouse Non-Rebreathing Circuits
- Includes 10 feet of 22mm EVAC hose and 1 f/air canister
- Can be used with standard activated charcoal canisters or passive waste gas evacuation systems

SCIVENA
SCIENTIFIC
ANESTHESIA INNOVATIONS

Formerly *LEI Medical*

(503) 822-1157

(877) 615-5632 toll-free

info@scivenascientific.com

12330 S.E. Highway 212
Clackamas, OR 97015