

# Slide Top Induction Chamber

P/N RES641


*Slide top design may help limit exposure from anesthetic gas exposure during subject removal*

**SCIVENA<sup>®</sup>**  
**SCIENTIFIC**  
ANESTHESIA INNOVATIONS



## Features:

- Ideal for mice, rats and other rodents
- Constructed of durable acrylic with anti-skid feet
- Precision machining slide top provides smooth opening
- Includes 15mm female inlet and 19/22mm male waste gas outlet



19/22MM waste gas  
evacuation port

15MM female  
gas inlet port

**Size: 1 Liter, 4" X 4" X 4" - For mice and other small rodents.**

**SCIVENA<sup>®</sup>**  
**SCIENTIFIC**

ANESTHESIA INNOVATIONS  
Formerly LEI Medical Since 1993

**(503) 655-2012**

(877) 615-5632 (toll Free)  
Info@ScivenaScientific.com  
12330 SE Hwy 212  
Clackamas, OR 97015